

**Bausch & Lomb**  
**Akreos® AO**

Micro Incision Lens

**INTRAOULAR LENSES**

DESCRIPTION: The Akreos Micro lens is a single-piece lens cut from a hydrophilic acrylic copolymer that contains UVA absorber. The overall length varies according to IOL power. The lens has a central optic zone and a peripheral zone which contains the UVA absorber.

INDICATIONS: The Akreos Micro posterior chamber lens is indicated for primary implantations for the visual correction of astigmatism in adult patients where a cataract has been removed or is not present. It is also indicated for secondary implantations after capsular bag replacement.


CONTRA-INDICATIONS: The Akreos Micro lens is contraindicated for patients with existing conditions, including the diagnosis or the treatment of a cataract, or previous eye surgery, unless the diagnosis or treatment has been completed.

WARNING: As with all surgical procedures, cataract surgery with the implantation presents risks of complications. These include, but are not limited to, infection (endophthalmitis), uveitis, macular edema, vitreous hemorrhage, retinal detachment, glaucoma, and hemorrhage. There is also a risk of damage to the lens capsule.

PRECAUTIONS FOR USE AND STORAGE: Store at room temperature. Avoid high temperatures (above 45°C) and direct sunlight. Do not expose the lens to heat or cold.

INDICATIONS: The Akreos Micro IOL is primarily packed in a pouch designed for the posterior chamber. The lens is supplied in a sterile pouch, which is held within the IOL. A patient card and self-adhesive labels are applied to provide individual identification information. The lens is supplied in its protective packaging and before opening the individual sterile pouch.

WARNING: All intraocular lenses are Kugelkernplastiken (gel-filled). They are not



can absorb substances that they contact (disinfectant, etc.). Consequently, Akreos lenses must be handled carefully to prevent damage to the lens. If the lens is damaged, it must be replaced, since the lens is bonded slightly behind the posterior capsule.

INDICATIONS: The Akreos Micro lens is indicated for primary implantations for the visual correction of astigmatism in adult patients where a cataract has been removed or is not present. The lens is supplied in a pouch designed for the posterior chamber.

CONTRA-INDICATIONS: The Akreos Micro lens is contraindicated for patients with existing conditions, including the diagnosis or the treatment of a cataract, or previous eye surgery, unless the diagnosis or treatment has been completed.

WARNING: All intraocular lenses are Kugelkernplastiken (gel-filled). They are not

designed to withstand temperatures above 45°C or direct sunlight. The lens is supplied in a pouch designed for the posterior chamber. The lens is supplied in a sterile pouch, which is held within the IOL. A patient card and self-adhesive labels are applied to provide individual identification information.

CONTRA-INDICATIONS:

WARNING:

AVERTISSEMENT:

AVIS:

INFORMACIÓN:

ADVARSEL:

VAROON:

AVVERTIMENTO:

AVERTIMENTO:

