
BAUSCH + LOMB

enVista ASPIRE™ HYDROPHOBIC ACRYLIC IOL

en
Devis Description
The enVista Aspire™ hydrophobic acrylic intraocular lens (IOL) (non-preloaded model: EA) was developed to replace the natural crystalline lens in adult patients in whom the cataractous lens has been removed. The composition and characteristics of the IOL are specified in the table below.

Indications

The enVista Aspire™ hydrophobic acrylic IOL (non-preloaded model EA) is indicated for primary implantation in the capsular bag of the eye in adult patients for the visual correction of aphakia following removal of a cataractous lens.

Warnings

- As with any surgical procedure, there is risk involved. Physicians considering IOL implantation under any of the following circumstances should weigh the potential risk/benefit ratio:
1. Recurrent severe anterior or posterior segment inflammation or uveitis.
2. Patients in whom the IOL may affect the ability to observe, diagnose, or treat posterior segment diseases.
3. Surgical difficulties at the time of cataract extraction, which might increase the potential for complications (e.g., persistent bleeding, significant iris damage, uncontrolled positive pressure, or significant vitreous prolapse or loss).
4. A distorted eye due to previous trauma or developmental defect in which appropriate support of the IOL is not possible.
5. Circumstances that would result in damage to the endothelium during implantation.
6. Suspected microbial infection.
7. Patients in whom neither the posterior capsule nor zonules are intact enough to provide support.

Physical Characteristics

Table with 2 columns: Lens/Haptic Material, Material Characteristics, Optic Type/Power, Dimensions, Spectral Transmittance. Details include hydrophobic acrylic, silicone copolymer, and various optical and material properties.

Precautions For Use And Storage

- 1. Do not attempt to sterilize the IOL as this can produce undesirable side effects.
2. Prior to opening, inspect vial pouch and vial for signs of damage that may affect integrity of device sterility. If damaged, do not use. The IOL should be used immediately after opening.
3. Do not use if product sterility or quality is thought to be compromised due to damaged packaging or signs of leakage (such as the loss of saline storage solution, or the presence of salt crystallization).
4. Store at room temperature. Do not freeze. Avoid high temperatures (>43°C / >109°F). Keep dry. Keep away from sunlight.
5. Do not soak or rinse the IOL with any solution other than sterile balanced salt solution or sterile normal saline.

- 6. Do not place the IOL in contact with surfaces where such contamination can occur.
7. Do not autoclave the IOL.
8. Do not re-use the IOL. It is intended for permanent implantation. If explanted, sterility and proper function cannot be ensured.
9. Neither the safety and effectiveness, nor the effects of the Aspire IOL optical design on depth of focus, contrast sensitivity, and subjective visual disturbances (glare, halo, etc.) have been evaluated clinically. MTF testing of the Aspire IOL optical design (used in model EA) (Figures 2 and 3) may aid the surgeon in understanding the theoretical image quality expected with the Aspire IOL compared to the enVista monofocal IOL MX60E. However, these do not fully assess all aspects of clinical difficulties under all conditions. Surgeons must weigh the potential benefits of the modified optical design of the Aspire IOL (model EA) against the potential for risks associated with a degradation in vision quality and the lack of clinical data to characterize the impact of the Aspire IOL optical design on contrast sensitivity and subjective visual disturbance. These considerations may be especially relevant to patients with certain pre-existing ocular conditions (prior corneal refractive surgery, irregular corneal astigmatism, severe corneal dystrophy, macular disease, optic nerve atrophy, etc.) or intraoperative conditions (posterior capsular rupture, complications in which the IOL stability could be compromised, inability to place IOL in capsular bag, etc).

10. The safety and effectiveness of the IOL have not been substantiated in patients with pre-existing ocular conditions and intraoperative complications (see below). Careful preoperative evaluation and sound clinical judgment should be used by the surgeon to decide the benefit/risk ratio before implanting an IOL in a patient with one or more of these conditions. Physicians considering IOL implantation in such patients should explore the use of alternative methods of aphakic correction and consider IOL implantation only if alternatives are deemed unsatisfactory in meeting the needs of the patient.

- 11. Patients with preoperative problems such as corneal endothelial disease, abnormal cornea, macular degeneration, retinal degeneration, glaucoma, and chronic drug miosis may not achieve the visual acuity of patients without such problems. The physician must determine the benefits to be derived from IOL implantation when such conditions exist.
12. A high level of surgical skill is required for IOL implantation. The surgeon should have observed and/or assisted in numerous implantations and successfully completed one or more courses on IOL implantation before attempting to implant the IOL.
13. As with any surgical procedure, there is risk involved. Potential complications accompanying cataract or limited surgery may include, but are not limited to the following: corneal endothelial damage, infection (endophthalmitis), retinal detachment, vitritis, cystoid macular edema, corneal edema, pupillary block, cystic membrane, iris prolapse, hypopyon, transient or persistent glaucoma, acute corneal decompensation, toxic anterior segment syndrome (TASS), and secondary surgical intervention. Secondary surgical interventions include, but are not limited to: IOL repositioning, lens replacement, vitreous aspiration or iridectomy for pupillary block, wound leak repair, and retinal detachment repair.
14. Care should be taken to remove viscoelastic from the eye at the close of surgery.

Before Surgery

- Retinal conditions or predisposition to retinal conditions, previous history of, or a predisposition to, retinal detachment or proliferative diabetic retinopathy, in which future treatment may be compromised by implanting this IOL.
• Amblyopia
• Clinically severe corneal dystrophy (e.g., Fuchs)
• Rubella, congenital, traumatic or complicated cataracts
• Extremely shallow anterior chamber, not due to swollen cataract

- Recurrent anterior or posterior segment inflammation of unknown etiology, or any disease producing an irritative reaction in the eye (e.g., iritis or uveitis)
• Aniridia
• Iris neovascularization
• Glaucoma (uncontrolled or controlled with medication)
• Microphthalmos or macrophthalmos
• Optic nerve atrophy
• Previous corneal transplant
• Pre-existing ocular conditions which may negatively impact stability of the implant

During Surgery

- Mechanical or surgical manipulation required to enlarge the pupil
• Vitreous loss (significant)
• Anterior chamber bleeding (significant)
• Uncontrollable positive intraocular pressure
• Complications in which the IOL stability could be compromised

Medical Device Re-Use Statement

If this product is reprocessed and/or re-used, Bausch + Lomb cannot guarantee the functionality, material structure, or cleanliness or sterility of the product. Re-use could lead to illness, infection and/or injury to the patient or user and, in extreme incidents, death. This product is labeled as 'single-use' which is defined as a device intended to be used once only for a single patient.

Calculation Of IOL Power Suggested A-Constant: 119.1 (OPTICAL BIOMETRY)

The recommended A-Constant is intended for use with axial length measurements obtained by optical biometry. Use of axial length measurements by other techniques (e.g., Aplanation A-scan) will normally require a different IOL constant. This number is a guideline only and is based on an evaluation of clinical data obtained using the IOL Master. The physician should determine preoperatively the power of the IOL to be implanted.

Directions For Use

- 1. Inspect vial pouch and vial for signs of damage that may affect integrity of device sterility. If damaged, do not use.
2. Prior to implanting, examine the IOL package for type, power, and proper configuration.
3. Open the peel pouch and remove the vial in a sterile environment.
4. Remove the lid from the vial.
5. Follow steps below.
a) With a pair of smooth forceps, remove the IOL from the vial by gently grasping the IOL haptic.
b) Rise the entire IOL with sterile balanced salt solution or sterile normal saline.
c) Examine the IOL thoroughly to ensure particles have not become attached to it, and examine the IOL optical surfaces for other defects.
d) The IOL may be soaked in sterile balanced salt solution until ready for implantation.
e) It is recommended to use an approved inserter per the Validated Inserters table below.

- 6. It is recommended to use an approved viscoelastic for lubrication of the IOL during implantation. See table below.
7. There are various surgical procedures that can be utilized, and the surgeon should select a procedure that is appropriate for the patient. Surgeons should verify that appropriate instrumentation is available prior to surgery.

Validated Inserters

Table with 3 columns: Model, Inserter, Viscoelastic. Lists models SA and SA and their corresponding inserters and viscoelastic materials.

Overview Of Clinical Studies

The enVista Aspire IOL model EA is a result of minor modifications (that did not require clinical study data) from the parent enVista model MX60. Clinical studies have not been conducted with the enVista Aspire IOL to assess the effect of its posterior aspheric surface on visual acuities, contrast sensitivity, or visual symptoms. Clinical studies have been conducted on the enVista single-piece IOL (model MX60) and the parent xact X-60 three-piece IOL (model X-60)*. The results of these studies are described herein.

* The AVS xact X-60 Intraocular Lens is not licensed for sale in Canada.

1. Summary Of Clinical Study For enVista Model MX60

A clinical study of the enVista Hydrophobic Acrylic Intraocular Lens, Model MX60, began in the United States on October 19, 2010. This prospective, single arm, open label study included a total of 122 subjects (122 eyes) at 6 clinical sites. Postoperatively, subjects underwent complete ophthalmic evaluations at regularly scheduled intervals through Form 4 (Postoperative Days 120-180).

Table 1 displays demographic information of subjects enrolled in the clinical trial. Table 2 displays BCVA results for best case subjects (those without clinically significant preoperative pathologies or macular degeneration at any time during the study) for 3 visits. At the Form 4 visit, 118 subjects (100%) achieved BCVA of 20/40 or better, which exceeds the FDA grid of 96.7%.

The key safety outcomes for this study are presented in Table 3. The rates of FDA defined potentially sight-threatening adverse events that occurred in the clinical trial at Form 4 were found to be less than the FDA Grid of Historical Controls. Two cumulative adverse events (2/122; 1.6%) of cystoid macular edema were reported through the Form 4 visit. One persistent adverse event (1/121; 0.8%) of cystoid macular edema was reported at the Form 4 visit. No serious ocular adverse events occurred during this study. One serious non-ocular adverse event of advanced leukemia with an outcome of death was reported during this study. The adverse event was determined by the study investigator to be unrelated to the investigational device, Model MX60 IOL. The incidence of adverse events experienced during the clinical trial was comparable to or lower than the incidence reported in the historic control ("FDA grid") population (see Table 6). The enVista IOL demonstrated favorable safety compared with the Control IOL and the historical control (cf. ISO 11979-7 SPE) population, with no increase in incidence or severity of adverse events (AEs) compared with the Control IOL and no serious adverse events (SAEs) in the study eye. Overall, no safety signals were associated with the IOLs during this study.

The results of clinical investigations provide reasonable assurance that the Model MX60 IOL is safe and effective for the visual correction of aphakia following cataract extraction.

Clinical Tables

Table 1: Subject Demographics. Table with 3 columns: Demographic, n, %. Lists gender, race, age, and other demographic data for the study population.

Table 2: Best Corrected Visual Acuity By Postoperative Visit (Best Case Analysis Set). Table with 4 columns: Visual Acuity, Form 2, Form 3, Form 4. Shows the percentage of subjects achieving various levels of visual acuity at different postoperative visits.

Table 3: ISO Defined Cumulative And Persistent Adverse Events Through Form 4 (Safety Analysis Set). Table with 5 columns: Adverse Event, n/N, %, ISO Grid, %, 99% CI, p-value. Lists various adverse events and their occurrence rates compared to ISO grid standards.

Table 4: Subject Demographics. Table with 2 columns: Demographic, n, %. Lists gender, race, and other demographic data for the study population.

Table 5: Cumulative Safety Events. Table with 5 columns: Adverse Event, n/N, %, ISO Grid, %, 99% CI, p-value. Lists cumulative safety events and their occurrence rates.

n: number of eyes reported with corresponding event. For cumulative event. For number of implanted eyes. For persistent event. No. number of eyes returned for the Form 4 examination with non-missing response for the corresponding adverse event. A subject could be reported with more than one AE.
Based on binomial distribution.
*Retained not for the null hypothesis H0: Present from study < Present from ISO Grid (per ISO 11979-7:2006 (E)).
**Missing at any time during the study.
†Present at Form 4.

Other Clinical Findings

All subjects in the safety analysis set were evaluated for IOL glistenings at Form 3 and Form 4 visits. IOL glistenings were evaluated via retroillumination slit lamp examination utilizing a photographic grading scale provided in the protocol. The grading scale consisted of (in order of severity): none, grade 0 (trace), grade 1, 2, 3, or 4. No glistenings of any grade were reported for any subject at any visit in the clinical study.

2. Summary Of Clinical Study For xact Model X-60 (Three-Piece IOL)

A clinical study of the xact Model X-60 IOL began in the United States on May 8th, 2002 and was conducted by Advanced Vision Science. A total of 383 subjects were enrolled, and 367 subjects were available for examination at one year, 312 were available at two years, and 281 were available at three years.

Table 4 displays demographic information of subjects enrolled in the clinical trial. Table 5 summarizes the best-corrected distance visual acuity (BCVA) results for best case subjects (those without clinically significant preoperative pathologies or macular degeneration at any time during the clinical trial).

Potentially sight threatening adverse events are listed in Table 6, along with the rate of occurrence in the clinical trial of the X-60 IOL, and are compared to the FDA Grid of Historical Controls. The number of patients included in the analysis of both cumulative and persistent adverse events in some cases was less than the number of patients who returned for examination and were available for analysis as a result of missing information in certain fields on the case report forms.

The results of clinical investigations provide reasonable assurance that the X-60 IOL is safe and effective for the visual correction of aphakia following cataract extraction.

Clinical Tables

Table 4: Subject Demographics. Table with 2 columns: Demographic, n, %. Lists gender, race, and other demographic data for the study population.

Table 5: Cumulative Safety Events. Table with 5 columns: Adverse Event, n/N, %, ISO Grid, %, 99% CI, p-value. Lists cumulative safety events and their occurrence rates.

n: number of eyes reported with corresponding event. For cumulative event. For number of implanted eyes. For persistent event. No. number of eyes returned for the Form 4 examination with non-missing response for the corresponding adverse event. A subject could be reported with more than one AE.
Based on binomial distribution.
*Retained not for the null hypothesis H0: Present from study < Present from ISO Grid (per ISO 11979-7:2006 (E)).
**Missing at any time during the study.
†Present at Form 4.

Table 5: Visual Acuity In Best Case Population

Table with 4 columns: Visual Acuity, 1 Year, 2 Years, 3 Years. Shows the percentage of subjects achieving various levels of visual acuity at different postoperative visits.

Table 6: Cumulative And Persistent Adverse Events

Table with 5 columns: Adverse Events, 1 Year, FDA Grid, 2 Years, 3 Years. Lists various adverse events and their occurrence rates compared to FDA grid standards.

IOL was exchanged due to patient complaint of blurred vision, despite good BCVA. Investigator suspected glistening might be related, however only modest improvement of vision was achieved after IOL exchange.
IOL was exchanged due to patient complaint of blurred vision, despite good BCVA. Investigator suspected glistening might be related, however only modest improvement of vision was achieved after IOL exchange.
IOL was exchanged due to patient complaint of blurred vision, despite good BCVA. Investigator suspected glistening might be related, however vision did not improve after IOL exchange. Since vision did improve after subsequent IOL exchange, the complaint of blurred vision was not associated with the IOL.

Other Clinical Findings

In the IDE clinical trial, "glistenings" were observed in some cases. Glistenings, known to sometimes occur in some other hydrophobic acrylic IOLs, are microscopic vacuoles within the optic of the IOL that are visible through the slit lamp as multiple small refractive specks. Analysis of the clinical data confirmed no effect of glistenings on visual outcomes.

Testing established that glistenings were eliminated by a change in the IOL hydration solution from 10.0% saline to 0.9% saline. This was confirmed in an additional clinical trial conducted outside of the United States. In this study, 172 eyes of 142 patients were examined at least once between 1 and 6 months, and 123 eyes of 101 patients were examined at least once between 6 months and 2 years. No glistenings were observed at any time.

Adverse Event Reporting

Adverse events and/or potentially sight-threatening complications that may be regarded as lens related and that were not previously expected in nature, severity or degree of incidence should be reported within five (5) days to Bausch & Lomb Incorporated. This information is being requested from all surgeons in order to document potential long-term effects of IOL implantation. Should any of these incidents take place in the European Union (EU), the event should also be reported to the competent authority of the EU member state in which the user is established.

Physicians are encouraged to report these events in order to aid in identifying emerging or potential problems with IOL. These problems may be related to a specific lot of lenses or may be indicative of long-term effects associated with these lenses or with IOLs in general. If you wish to report a problem in the USA, please call Bausch + Lomb at 1-800-338-2020. Outside the USA, contact information can be found on www.bausch.com/contact. In case of device explant/extraction from a patient, keep the device and contact Customer Service for return instructions.

How Supplied

The enVista IOL (model EA) is individually packaged in a sterile vial (containing blank bank saline) within a peel pouch, and should only be opened under sterile conditions.

A patient card and self-adhesive labels are supplied to provide traceability of the lens. The package is sterilized by gamma irradiation.

Expiry Date

Sterility is guaranteed unless the pouch is damaged or opened. The expiry date on the lens package is the sterility expiry date. This lens should not be implanted after the indicated sterility expiry date.

Safe Disposal

Dispose/discard of the unused or contaminated device/equipment and/or packaging by following applicable safe disposal procedures, and in accordance with applicable laws and regulations regarding the disposal of biohazardous materials.

Patient Registration Instructions And Reporting Registration

Each patient who receives an enVista IOL must be registered with Bausch + Lomb at the time of lens implantation. Registration is accomplished by completing the Implant Registration Card that is enclosed in the lens package and mailing it to Bausch + Lomb. Patient registration is essential and will assist Bausch + Lomb in responding to adverse reaction reports and/or potentially sight-threatening complications. An implant identification card is supplied in the lens package and must be given to the patient.

Symbols And Abbreviations Used On Labeling

Table with 2 columns: Symbol or Abbreviation, Symbol or Abbreviation Title. Lists various symbols and their corresponding titles used on the product labeling.

Bausch + Lomb Incorporated
1400 North Goodman Street
Rochester, NY 14609 USA

Manufactured by: Bausch + Lomb Incorporated
21 Park Place Blvd, North Clearwater, FL 33759 USA
®/TM are trademarks of Bausch + Lomb Incorporated or its affiliates.
© 2023 Bausch + Lomb Incorporated or its affiliates
www.bausch.com/symbols
Rev. 2023-09
4172100

Each patient who receives an enVista IOL must be registered with Bausch + Lomb at the time of lens implantation. Registration is accomplished by completing the Implant Registration Card that is enclosed in the lens package and mailing it to Bausch + Lomb. Patient registration is essential and will assist Bausch + Lomb in responding to adverse reaction reports and/or potentially sight-threatening complications. An implant identification card is supplied in the lens package and must be given to the patient.

Figure 2



Figure 3

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