Stellaris Elite™

vision enhancement system

REF BL15455

en

Operator's Manual Addendum

THIS ADDENDUM IS INTENDED TO PROVIDE USER INFORMATION RELATED TO THE PHOTOCOAGULATION LASER FOR POSTERIOR SEGMENT SURGERY. THE OPERATOR SHOULD THOROUGHLY READ AND UNDERSTAND ALL SAFETY PRECAUTIONS AND WARNINGS IN THE OPERATOR'S MANUAL PRIOR TO USING THE SYSTEM.

THE *Stellaris Elite*™ vision enhancement system WITH SOFTWARE VERSION 5.25 AND ABOVE LIMITS THE LASER MAXIMUM POWER OUTPUT SETTING TO 1200 mW FOR EVERY REGION WHERE THE 2000 mW LASER IS REGISTERED, THIS SOFTWARE OPTION WILL BE AUTOMATICALLY ENABLED TO EXTEND THE MAXIMUM POWER OUTPUT SETTING TO 2000 mW. DISABLING THE LASER POWER EXTENSION SOFTWARE OPTION LIMITS THE *Stellaris Elite*™ vision enhancement system TO OPERATE AT A MAXIMUM LASER POWER OF 1200 mW.

Changes to the User Interface:

System configuration screen:

- An Extended Maximum Laser Power software option is added to the system configuration screen. (See *Figure 1*).
- In all systems where a 2000 mW laser is registered, this software option will automatically be enabled upon software update.
- If the software option is disabled at any time, contact your service engineer to re-enable the software option.
- For all new systems, the software option is factory installed before shipment where a 2000 mW laser is registered (not available for systems built for China).

Note: Performing system upgrade to version 5.25 and above would set the default laser limit to 1200 mW. To extend the maximum power output setting to 2000 mW, the laser power extension software option must be enabled. If 2000 mW is the currently registered limit in your region, this will automatically be enabled. Your service engineer will be able to enable this option.

Figure 1: System Configuration Screen with Laser Power Extension software option

www.bausch.com/ifu

Bausch & Lomb Incorporated
1400 North Goodman Street
Rochester, NY 14609 USA

Bausch & Lomb GmbH
Brunsbütteler Damm 165/173
13581 Berlin, Germany

UK Responsible Person
Bausch & Lomb U.K. Ltd.
106 London Road
Kingston upon Thames
KT2 6TN UK

Manufactured by:
Bausch & Lomb Incorporated
3365 Tree Court Industrial Blvd.
St. Louis, MO 63122 USA

www.bausch.com/symbols

STELLARIS ELITE is a trademark of Bausch & Lomb Incorporated or its affiliates.

© 2023 Bausch & Lomb Incorporated or its affiliates

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