BAUSCH + LOMB

enVista

HYDROPHOBIC ACRYLIC IOL

ASPIRE™

en

Device Description

The enVista Aspire™ hydrophobic acrylic intraocular lens (IOL) (non-preloaded model: EA / preloaded models: EPA, EC20A, EC24A) was developed to replace the natural crystalline lens in adult patients in whom the cataractous lens has been removed. The composition and characteristics of the IOL are specified in the table below.

All IOL optical designs are associated with a certain amount of depth of focus. The Aspire IOL is classified as a monofocal IOL. Monofocal IOLs provide a limited depth of focus. The Aspire IOL uses an optical modification of the posterior aspheric surface to create a small continuous increase in IOL power within the central 1.5 mm diameter to slightly extend the depth of focus. However, clinically meaningful extension of the depth of focus has not been demonstrated in clinical trials.

In general, extending the depth of focus negatively affects the quality of vision at far distances. For the Aspire IOL, the power profile decreases towards the periphery outside the central 1.5 mm optic diameter, resulting in comparable distance image quality to the enVista monofocal IOL MX60E for a 3 mm pupil. Image quality has been estimated using non-clinical bench testing in an ISO2:±0.28 µm SA cornea model eye that simulates the spherical aberration from the cornea (see Figures 2 and 3).

The optic is designed with the SureEdge™ posterior-squared step edge to provide a 360-degree PCO barrier. The IOL employs an Accuset™ haptic with a broad, modified C-loop design and optic-haptic offset to facilitate improved contact and stability within the capsular bag. The enVista material that makes up the TruSight™ optic has been assessed for glistening-free capacity and scratch resistance. Incorporated StableFlex™ technology allows for enhanced IOL compliance for ease of loading, premium control through IOL delivery, and efficient optical recovery following implantation.

Physical Characteristics

Lens / Haptic Material	Hydrophobic acrylic (hydroxyethyl methacrylate (HEMA)-polyethylene glycol phenyl ether acrylate (poly(EG)PEA)-styrene copolymer, crosslinked with ethylene glycol dimethacrylate)
Material Characteristics	Index Of Refraction: 1.53 @ 35°C; Specific Gravity: 1.19 g/ml
Optic Type / Powers	Aspheric / +6.0 to +34.0 Diopters in 0.5 Diopter increments (SE – Spherical Equivalent)
Dimensions	Body Diameter: 6.0 mm; Overall Diameter: 12.5 mm; Haptic Angle: 0°
Spectral Transmittance	Ultraviolet: UV (389) 10% transmittance for +20.0 Diopter IOL. See Figure 1 with chart's X value = Wavelength (nm) and Y value = % Transmittance; chart compares the transmittance curve of an enVista IOL to a 53-Year-Old Human Lens. A = +20 Diopter enVista EE IOL and B = 53-Year-Old Human Lens. NOTE: Light transmittance values for an IOL material may vary slightly depending on the method of measurement. Reference: 53-Year-Old Human Lens data from Boettner, E.A. and Wolter, J.R., "Transmission of the Ocular Media," Investigative Ophthalmology, 1:776-783, 1962.
MTF Curves	See Figure 2 with chart's X value = Focus Position (D) and Y value = MTF at 50 lp/mm; chart compares the through focus MTF at 50 lp/mm, 546 nm light in ISO2 model eye with 3 mm aperture A = enVista Aspire IOL and B = +20 Diopter enVista MX60E IOL See Figure 3 with chart's X value = Frequency (cycles/mm) and Y value = MTF; chart compares the through frequency MTF 546 nm light in ISO2 model eye with 3 mm aperture A = enVista Aspire IOL and B = +20 Diopter enVista MX60E IOL

Indications

The enVista Aspire hydrophobic acrylic IOL (non-preloaded model: EA / preloaded models: EPA, EC20A, EC24A) is indicated for primary implantation in the capsular bag of the eye in adult patients for the visual correction of aphakia following the removal of a cataractous lens for improved uncorrected distance vision.

enVista Aspire IOLs are intended to be folded and placed in the posterior chamber of the eye after the cataractous crystalline lens has been removed.

Clinical Benefits

The enVista Aspire posterior chamber IOL provides a safe and effective product for the visual correction of aphakia secondary to the removal of the cataractous lens in adult patients.

Contra-Indications

Implantation is not advisable when the IOL may aggravate an existing condition, interfere with the diagnosis or the treatment of a pathology, or present a risk to the sight of the patient. These conditions include, but are not limited to: uncontrolled glaucoma, rubeotic cataract, retinal detachment, atrophy of the iris, microphthalmia, developing chronic eye infections, endothelial corneal dystrophy, peri-operative complications (such as vitreous loss, hemorrhage), or foreseeable postoperative complications. Additional considerations are listed in the Warnings and Precautions for Use and Storage sections.

Warnings

As with any surgical procedure, there is risk involved. Physicians considering IOL implantation under any of

the following circumstances should weigh the potential risk/benefit ratio:

- Recurrent severe anterior or posterior segment inflammation or uveitis.
- Patients in whom the IOL may affect the ability to observe, diagnose, or treat posterior segment diseases.
- Surgical difficulties at the time of cataract extraction, which might increase the potential for complications (e.g., persistent bleeding, significant iris damage, uncontrolled positive pressure, or significant vitreous prolapse or loss).
- A distorted eye due to previous trauma or developmental defect in which appropriate support of the IOL is not possible.
- Circumstances that would result in damage to the endothelium during implantation.
- Suspected microbial infection.
- Patients in whom neither the posterior capsule nor zonules are intact enough to provide support.

Precautions For Use And Storage

- Do not attempt to resterilize the IOL as this can produce undesirable side effects.
- Prior to opening, inspect vial pouch and vial for signs of damage that may affect integrity of device sterility. If damaged, do not use. The IOL should be used immediately after opening.
- Do not use if product sterility or quality is thought to be compromised due to damaged packaging or signs of leakage (such as the loss of saline storage solution, or the presence of salt crystallization).
- Store at room temperature. Do not freeze. Avoid high temperatures (>43°C / >109°F). Keep dry. Keep away from sunlight. Do not use if the packaging is exposed to environmental conditions outside of those specified.
- Do not soak or rinse the IOL with any solution other than sterile balanced salt solution or sterile normal saline.
- Do not place the IOL in contact with surfaces where such contamination can occur.
- Do not autoclave the IOL.
- Do not re-use the IOL. It is intended for permanent implantation. If explanted, sterility and proper function cannot be ensured.
- Neither the safety and effectiveness, nor the effects of the enVista Aspire IOL optical design on depth of focus, contrast sensitivity, and subjective visual disturbances (glare, halo, etc.) have been evaluated clinically.
- The safety and effectiveness of the IOL have not been substantiated in patients with pre-existing ocular conditions and intraoperative complications (see below). Careful preoperative evaluation and sound clinical judgment should be used by the surgeon to decide the benefit/risk ratio before implanting an IOL in a patient with one or more of these conditions. Physicians considering IOL implantation in such patients should explore the use of alternative methods of aphakic correction and consider IOL implantation only if alternatives are deemed unsatisfactory in meeting the needs of the patient.

Before Surgery

- Retinal conditions, or predisposition to retinal conditions, previous history of, or a predisposition to retinal detachment, or proliferative diabetic retinopathy in which future treatment may be compromised by implanting this IOL
- Amblyopia
- Clinically severe corneal dystrophy (e.g., Fuchs')
- Rubella, congenital, traumatic or complicated cataracts
- Extremely shallow anterior chamber, not due to swollen cataract
- Recurrent anterior or posterior segment inflammation of unknown etiology, or any disease producing an inflammatory reaction in the eye (e.g., iritis or uveitis)
- Aniridia
- Iris neovascularization
- Glaucoma (uncontrolled or controlled with medication)
- Microphthalmos or macrophthalmos
- Optic nerve atrophy
- Previous corneal transplant
- Pre-existing ocular conditions which may negatively impact stability of the implant

During Surgery

- Mechanical or surgical manipulation required to enlarge the pupil
 - Vitreous loss (significant)
 - Anterior chamber bleeding (significant)
 - Uncontrollable positive intraocular pressure
 - Complications in which the IOL stability could be compromised
- Patients with preoperative problems, such as corneal endothelial disease, abnormal cornea, macular

- degeneration, retinal degeneration, glaucoma, and chronic drug miosis may not achieve the visual acuity of patients without such problems. The physician must determine the benefits to be derived from IOL implantation when such conditions exist.
- A high level of surgical skill is required for IOL implantation. The surgeon should have observed and/or assisted in numerous implantations and successfully completed one or more courses on IOL implantation before attempting to implant the IOL.
 - As with any surgical procedure, there is risk involved. Potential complications accompanying cataract or implant surgery may include, but are not limited to, the following: corneal endothelial damage, infection (endophthalmitis), retinal detachment, vitritis, cystoid macular edema, corneal edema, pupillary block, cyclitic membrane, iris prolapse, hypopyon, transient or persistent glaucoma, acute corneal decompensation, toxic anterior segment syndrome (TASS), and secondary surgical intervention. Secondary surgical interventions include, but are not limited to: IOL repositioning, IOL replacement, vitreous aspiration or iridectomy for pupillary block, wound leak repair, and retinal detachment repair.
 - Care should be taken to remove all viscoelastic from the eye at the close of surgery.

After Surgery

- The patient should meet with the surgeon soon after surgery to confirm proper placement of the IOL.

Medical Device Re-Use Statement

If this product is reprocessed and/or re-used, Bausch + Lomb cannot guarantee the functionality, material structure, or cleanliness or sterility of the product. Re-use could lead to illness, infection and/or injury to the patient or user and, in extreme incidents, death. This product is labeled as 'single-use' which is defined as a device intended to be used once only for a single patient.

Calculation Of Lens Power

Suggested A-Constant: 119.1 (OPTICAL BIOMETRY)

The recommended A-Constant is intended for use with axial length measurements obtained by optical biometry. Use of axial length measurements by other techniques (e.g., Applanation A-scan) will normally require a different lens constant. This number is a guideline only and is based on an evaluation of clinical data obtained using the IOL Master. The physician should determine preoperatively the power of the lens to be implanted.

Directions For Use

- Inspect vial pouch and vial for signs of damage that may affect integrity of device sterility. If damaged, do not use.
- Prior to implanting, examine the IOL package for type, power, and proper configuration.
- Open the peel pouch and remove the vial in a sterile environment.
- Remove the lid from the vial.
- Follow steps below.
 - Non-preloaded IOL (model: EA)
 - With a pair of smooth forceps, remove the IOL from the vial by gently grasping the IOL haptic.
 - Rinse the entire IOL with sterile balanced salt solution or sterile normal saline.
 - Examine the IOL thoroughly to ensure particles have not become attached to it, and examine the IOL optical surfaces for other defects.
 - The IOL may be soaked in sterile balanced salt solution until ready for implantation.
 - It is recommended to use an approved inserter per the Validated Inserters table below.
 - Preloaded lens (models: EPA, EC20A, EC24A)
 - Remove the SnapSet™ IOL shuttle from the vial by grasping the tab on the top.
 - The SnapSet IOL shuttle may be flushed with sterile balanced salt solution during preparation for delivery.
 - The SnapSet IOL shuttle preloaded with the enVista IOL (models: EPA, EC20A, EC24A) is designed for use with the EyeGility™ inserter for enVista preloaded.
- It is recommended to use an approved viscoelastic for lubrication of the IOL during implantation. See table below.
- There are various surgical procedures that can be utilized, and the surgeon should select a procedure that is appropriate for the patient. Surgeons should verify that appropriate instrumentation is available prior to surgery.
- The IOL and insertion device should be discarded if the IOL has been held in the folded state within the insertion device for more than 20 minutes. Not doing so may result in damage to the IOL.

Validated Inserters

Model	inserter	Viscoelastic
EA	BLIS (BLIS-X1 cartridge)	Amvisc™ Plus, Amvisc™, OcuCoat™
	INJ100	
EPA (IOL only)	EyeGility (2.0 mm or 2.4 mm)	
EC20A (IOL + inserter)	EyeGility 2.0 mm	
EC24A (IOL + inserter)	EyeGility 2.4 mm	

Overview Of Clinical Studies

Clinical studies have been conducted on the enVista single-piece IOL (model: MX60) and the parent xact X-60 three-piece IOL (model X-60). The results of the clinical study provide reasonable assurance that the model MX60 IOL is safe and effective for the visual correction of aphakia following cataract extraction.

The enVista hydrophobic acrylic intraocular lens, model MX60, is the parent lens of the models listed in the device description. For the enVista MX60 clinical study, all subjects in the safety analysis set were evaluated for IOL glistenings at Form 3 and Form 4 visits. IOL glistenings were evaluated via retroillumination slit lamp examination utilizing a photographic grading scale provided in the protocol. The grading scale consisted of (in order of severity), "none, grade 0 (trace), grade 1, 2, 3, or 4." No glistenings of any grade were reported for any subject at any visit in the clinical study.

Adverse Event Reporting

Adverse events and/or potentially sight-threatening complications that may be regarded as lens related and that were not previously expected in nature, severity, or degree of incidence should be reported within five (5) days to Bausch + Lomb Incorporated. This information is being requested from all surgeons in order to document potential long-term effects of IOL implantation. Should any of these incidents take place in the European Union (EU), the event should also be reported to the competent authority of the EU member state in which the user is established.

Physicians are encouraged to report these events in order to aid in identifying emerging or potential problems with the IOL. These problems may be related to a specific lot of lenses or may be indicative of long-term effects associated with these lenses or with IOLs in general. If you wish to report a problem, in the USA, please call Bausch + Lomb at 1-800-338-2020. Outside the USA, contact information can be found on www.bausch.com/contactus. In case of device explant/extraction from a patient, keep the device and contact Customer Service for return instructions.

Summary Of Safety And Clinical Performance (SSCP)

The Basic UDI-DI is the key identifier of a device in the EUDAMED (European Database on Medical Devices), on certificates and on the Declaration of Conformity included in the SSCP provided by Bausch + Lomb Incorporated and can be accessed by following this link: <https://ec.europa.eu/tools/eudamed>. The role of the UDI-DI is to unambiguously identify a device in the distribution chain; it's the primary identifier of a device model and is composed of a unique numeric or alphanumeric code.

How Supplied

Non-preloaded enVista IOL:
The enVista Aspire IOL (model: EA) is individually packaged in a sterile vial (containing blood bank saline) within a peel pouch, and should only be opened under sterile conditions.

Preloaded IOL:
The enVista Aspire IOL (models: EPA, EC20A, EC24A) is preloaded in the SnapSet IOL shuttle and individually packaged in a sterile vial (containing blood bank saline), within a peel pouch, and should only be opened under sterile conditions.

A patient card and self-adhesive labels are supplied to provide traceability of the lens. The package is sterilized by gamma irradiation.

Expiry Date

Sterility is guaranteed unless the pouch is damaged or opened. The expiry date on the lens package is the sterility expiry date. This lens should not be implanted after the indicated sterility expiry date.

Safe Disposal

Dispose/discard of the unused or contaminated device/ equipment and/or packaging by following applicable safe disposal procedures, and in accordance with applicable laws and regulations regarding the disposal of biohazardous materials.

Patient Registration Instructions And Reporting Registration

Each patient who receives an enVista IOL must be registered with Bausch + Lomb at the time of lens implantation. Registration is accomplished by completing the Implant Registration Card that is enclosed in the lens package and mailing it to Bausch + Lomb. Patient registration is essential and will assist Bausch + Lomb in responding to adverse reaction reports and/or potentially sight-threatening complications. An implant identification card is supplied in the lens package and must be given to the patient.

Symbols And Abbreviations Used On Labeling

Symbol or Abbreviation	Symbol or Abbreviation Title
SE	Spherical Equivalent
CYL	Cylinder
IOL	Intraocular Lens
PC	Posterior Chamber
UV	Ultraviolet
D	Diopter
Ø _B	Body Diameter (Optic Diameter)
Ø _T	Overall Diameter (Overall Length)
R _X ONLY	Caution: Federal (US) law restricts this device to sale by or on the order of a physician
	Fee paid for waste management
	Monofocal Intraocular Lens
	Preloaded Monofocal Intraocular Lens

Figure 1


Figure 2



Figure 3

Bausch + Lomb Incorporated
1400 North Goodman Street
Rochester, NY 14609 USA

Bausch + Lomb GmbH
Brunsbütteler Damm 165/173
13581 Berlin, Germany

2797

Manufactured by:
Bausch + Lomb Incorporated
21 Park Place Blvd. North
Clearwater, FL 33759 USA

Bausch + Lomb Netherlands B.V.
Koolhovenlaan 110, 1119 NH Schiphol-Rijk
The Netherlands

Bausch + Lomb Swiss AG
Industriestrasse 15A
6301 Zug, Switzerland

®/TM are trademarks of
Bausch + Lomb Incorporated
or its affiliates.
© 2024 Bausch + Lomb Incorporated
or its affiliates

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Bausch & Lomb GmbH
Bausch & Lomb Swiss AG